



Announcing a
NEW
BPA
Agreement

EuroAmerican IP, LLC
Specialized Medical Products



Solution Based Distribution

GRIP-LOK[®]

Safe. Simple. Secure.



- **BPA Number:**
VA199-17-A-0029
- **FSS/VA Contract Number:**
V797P-2181D
- **Effective Dates:**
November 1, 2016 - October 31, 2017
- **Item Number:**
GL-3400LFC-BX
- **Description:**
Grip-Lok® Universal
Large Foley Catheter
Securement Device
- **Quantity Per Box:**
100 pieces
- **Price Breaks:**

| QUANTITY (Boxes) | COST Per Box |
|-----------------------------|-------------------------|
| 1-500 | \$290.04 |
| 501-1000 | \$287.39 |
| 1001-1500 | \$283.57 |
| 1501+ | \$279.67 |

CONTACT

EuroAmerican IP, LLC **Specialized Medical Products**



Solution Based Distribution

EuroAmerican IP, LLC
2511 NW 16th Lane, Bay #2
Pompano Beach, Florida 33064

Toll Free: 866-972-6467
Local: 954-977-2929

Info@EuroAmericanProducts.com
www.EuroAmericanProducts.com

Debra J. Paul, Owner / CEO
Debra@EuroAmericanProducts.com

Tami Dana, Co-Owner / COO
Tami@EuroAmericanProducts.com



Application of **GL-3400LFC-BX** for Large Foley Catheter Securement

Prepare the skin according to the standard hospital protocol for dressing application. Skin prep or hair removal may be required on some patients for better adhesion. Skin must be **CLEAN** and **DRY** prior to Grip-Lok® application.



1

Open the top fabric over-strap and position the catheter between the foam strip and the strap so that the catheter bifurcation crotch is centered. Pull the paper backing from one side of the Grip-Lok® and secure to skin, then the other.



2

Separate the two sides of the fabric over-strap.



3

Secure one side of the top strap over the catheter at the bifurcation crotch.



4

Secure the other side of the top strap over the drainage lumen. The catheter is secured and the degree of rotational, and limited horizontal movement can be adjusted by applying the top fabric over-strap in a tighter fashion.



GL-3400LFC-BX will also fit triple lumen catheters.

For more information
visit grip-lok.com

